

RESMED



Join Wake Up to Sleep—ResMed's sleep apnea support community—at [WakeUpToSleep/Join](#).


 **ResShop**
better sleep, better life.
Click Here To Buy

Fitting

Swift FX Bella Series



1 Hold the pillows at your nose, ensuring the ResMed logo is at the top. Pull one strap around and behind your ear.



2 Ensure that the strap lies flat and comfortably across the ear. Take care to avoid twisting the strap.



3 Repeat the same process with the other strap.



4 Adjust the mask's fit by pulling both of the upper straps at the same time, along the notches, until you get a comfortable fit.



5 Repeat the process on the mask's lower straps. Do not overtighten.



6 Check that the pillows are upright at the entrance to the nose, then connect the free end of the short tube assembly to the air tubing of your device. Now it's ready for use.

OPTIONAL



You also have the option to wear your Swift FX Bella nasal pillows with the included Swift FX or Swift FX for Her headgear. For more information visit [ResMed.com](#).



Combine a Swift FX Bella Series mask with ResMed's S9 AutoSet™ therapy device for a comfortable, quiet and lifestyle-friendly system.



Swift™ FX Bella Series
NASAL PILLOWS SYSTEM

Choice fit. Complete comfort.

RESMED

Choice fit. Complete comfort.



Swift™ FX Bella Series
NASAL PILLOWS SYSTEM

[ResMed.com](#)

Featuring a minimal design that provides a feeling of freedom, the **Swift™ FX Bella Series** is lighter, softer and more streamlined for maximum comfort.

*Combined responses of those expressing a preference in an independent patient survey in the USA, UK, Germany, France & Australia, from September 2012. Visit [www.resmed.com/maskbrand](#).

ResMed Corp San Diego, CA, USA +1 858 836 5000 or 1 800 424 0737 (toll free). ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). See [ResMed.com](#) for other ResMed locations worldwide. AutoSet, S9, S9 AutoSet and Swift are trademarks of ResMed Ltd and are registered in the U.S. Patent and Trademark Office. ©2013 ResMed. Specifications may change without notice. 1015999/1 2013-01

Swift FX Bella Series


Easy, comfortable style

With its unobtrusive design and sophisticated cushion technology, Swift FX Bella Series, nasal pillows provide a comfortable fit with a feeling of freedom for men and women.

Pillows (3 sizes) come packaged with the mask system to fit a wide range of users.



Swift FX Bella Series



Bella loops
fit comfortably around the ears and stay out of the hair for a more streamlined look and feel

Dual-wall nasal pillows
provide a gentle, comfortable seal on contact

Soft, narrow cushion with integrated flexible chamber
makes side sleeping comfortable and helps maintain seal during movement

Ultra smooth 360° rotating elbow
rotates with easy, fluid movement

Quiet vents
diffuse and direct airflow quietly and gently away from the user and bed partner


Spring-flex lightweight tube
stretches and bends easily to reduce pull on the mask and enhance stability

"Lock and key" feature
allows for easy attachment and detachment of Bella loops (or Swift FX headgear straps)

Size and name indicators
make it easy to identify the size as well as the correct orientation of the mask

Flexible pillows with trampoline action
accommodate movement in multiple directions while maintaining seal. Three pillow sizes come packaged with mask to fit most patients

Swift FX Bella components



Pillows
61520 XS
61521 S
61522 M
61523 L¹

Bella loops accessory
61581 (1 pair)

Swift FX for Her headgear assembly
(headgear, pink backstrap and pink soft wraps)
61543

Frame assembly
(short tube, elbow and swivel)
61528

Frame system
(no headgear)
61510 XS
61511 S
61512 M
61513 L¹

Swift FX for Her headgear and Bella loops accessory combo pack
61580

PRODUCT CODES

Complete Systems

Swift FX Bella Nasal Pillows Mask

Style	Part Number
Swift FX Bella	61560


Swift FX Bella complete system includes Bella loops attached to the pillows frame, Swift FX for Her headgear with pink backstrap and pink soft wraps. Also included are Extra Small (XS), Small (S) and Medium (M) pillows.¹

Reimbursement Codes

HCPCS	Code Descriptor
A7034	Mask interface—mask or cannula, 1 per 3 months
A7033	Pillow for use on nasal cannula type interface—replacement, 2 pairs per month
A7035	Headgear used with positive airway pressure device, 1 per 6 months

¹ Large pillows are available separately.

Swift FX Bella Gray components



Pillows
61520 XS²
61521 S
61522 M
61523 L

Bella loops accessory
61582 (1 pair)

Swift FX headgear assembly
(headgear, blue backstrap and gray soft wraps)
61529

Frame assembly
(short tube, elbow and swivel)
61528

Frame system
(no headgear)
61510 XS²
61511 S
61512 M
61513 L

Swift FX headgear and Bella Gray loops accessory combo pack
61583

PRODUCT CODES

Complete Systems

Swift FX Bella Gray Nasal Pillows Mask

Style	Part Number
Swift FX Bella Gray	61568

Swift FX Bella Gray complete system includes Bella Gray loops attached to the pillows frame, Swift FX headgear with blue backstrap and gray soft wraps. Also included are Small (S), Medium (M) and Large (L) pillows.²

Reimbursement Codes

HCPCS	Code Descriptor
A7034	Mask interface—mask or cannula, 1 per 3 months
A7033	Pillow for use on nasal cannula type interface—replacement, 2 pairs per month
A7035	Headgear used with positive airway pressure device, 1 per 6 months

² Extra Small (XS) pillows are available separately.