

XT Heated Humidifier

Instruction Manual

APEX MEDICAL S.L.

Alameda de Recalde 50, Despacho 304
48008 Bilbao, Vizcaya, SPAIN

9, Min Sheng St., Tu-Cheng,
Taipei County, 236, Taiwan

USER'S MANUAL

English

User's manual P. 1

Español

Manual de Instrucciones P.14

Português

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Français

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German

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German

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IMPORTANT SAFEGUARDS – PLEASE MAKE THIS PRECAUTIONS LIST AVAILABLE TO ALL USERS and PROVIDERS

READ ALL INSTRUCTIONS BEFORE USING THIS SYSTEM

DANGER – Reducing the risk of electrocution:

- Always unplug this product immediately after using.
- Do not use while bathing.
- Do not place or store product where it can fall into or be pulled into a tub or sink.
- Do not place in or immerse into water or other liquids.
- Do not reach for the unit if it is immersed in water. Unplug the unit immediately.

WARNING – Reducing the risk of burns, electrocution, fire, or personal injury.

- Do not operate CPAP (Continuous Positive Airway Pressure, which is used for Obstructive Sleep Apnea treatment) before mask and humidifier is installed.
- Disconnect the humidifier from CPAP unit when not in use. Water entering the CPAP unit may result in electrical hazard (shock) or damage to the CPAP unit.
- Do not use if water comes in contact with the CPAP unit or enters the tubing.
- XT Heated Humidifier (9S-006) is not for use with patients whose supraglottic airways have been bypassed.
- Do not use the XT Heated Humidifier (9S-006) on the floor.
- The docking platform of the XT Heated Humidifier is only intended for the 9S-005 (XT) series CPAP. If used with other devices, you may encounter potential electrical hazards (shock) or damage to the XT Heated Humidifier unit.
- Do not cover the humidifier with blankets, clothes or other material.
- Do not remove the water chamber while the humidifier is on and heating the chamber. Do not touch the heating plate for 30 minutes after the device is disconnected from the electrical socket.
- This product should never be left unattended when plugged in.
- Close supervision is necessary when this product is used by, on or near children or individuals.
- Use this product only for its intended use as described in this manual. Use this product only under a physician's supervision. Do not use the equipment with accessories other than those supplied with this unit. Contact your local dealer/EU representative for additional information.
- Never operate this product if a) it has a damaged cord or plug, b) it is not operating properly, c) it has been dropped or otherwise traumatized, d) it has been immersed into water. Return the product to a qualified service center, distributor or EU representative for examination and repair.
- Keep the cord away from heated surfaces.
- Never block the air openings of this product or allow objects to fall or be inserted into the air vent openings. Do not place unit on soft surfaces such as beds or couches, where the air openings may become obstructed.
- This device IS NOT to be used if the room temperature is warmer than 35°C (95°F) to prevent the temperature of air delivered to nasal mask from exceeding 41°C (104°F).

CAUTION –

- Use humidifier only with an XT series CPAP system and CE certified or 510(k) certified accessories

NOTE –

- U.S. Federal law restricts this device to sale by, or on the order of, a licensed physician.

NOTES, PRECAUTIONS AND WARNINGS

- DANGER:** Indicates an extreme and immediate hazardous situation which can result in death or serious injury.
- WARNING:** Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury
- CAUTION:** Refers to correct operating or maintenance procedures that is required in order to prevent damage or destruction to the equipment or other property.
- NOTE:** Refers to information that requires special attention.

SYMBOLS

Manufacturer

Authorised representative in the European community

The BF symbol indicates that this product complies with regulations protecting against electrical shock for type BF equipment.

Attention to warnings or safety precautions.

Grounding terminal

High Water Mark Fill line (for water chamber)

Air flow direction

Disposal of Electrical & Electronic Equipment (WEEE):

This product should be handed over to an applicable collection point for the recycling of electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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1. Introduction

This manual is provided to explain initial set up and operating instructions. This manual should be stored in a safe location and be readily accessible for future reference.

1.1 General Information

This system has been tested and successfully approved by the following standards.

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices for IEC 60601-1-2:2004. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This unit generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the unit.
- Increase the distance between the conflicting equipment.
- Connect the unit into an outlet that is on a different circuit than the other device(s).
- Consult the manufacturer or field service technician for help.


This system has been tested and is in compliance with the following volunteer standards.
FDA

1.2 Intended use

XT Heated Humidifier is intended for use with the XT series CPAP unit. It has been designed to increase the moisture in the air delivered from the CPAP, alleviating symptoms of a dry nose and throat which may result from ordinary airflow from the CPAP unit.

NOTE: This equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen and/or nitrous oxide.

2. Product description

- Water Chamber:** The water chamber contains the water that will be humidified into the airflow. The base plate can be removed to clean the chamber.
- Fill Line:** This symbol is located on the clear plastic of the chamber and is the high water mark that indicates the maximum amount of water you can put into the chamber.
- Flow Direction Indicator:** This is an arrow symbol on top of the inlet port that indicates the air flow direction
- Baffle:** The baffle is used to prevent sloshing and water spillage when the unit is moved or inadvertently bumped.
- Extraction Tool:** It is used to separate the base of the water chamber from the chamber itself (for cleaning).
- Heater Plate:** Heats the water in the water chamber.
- Locking Tab:** Keeps the CPAP in position on the humidifier platform.
- Heater Plate Switch:** When the chamber is in place, it will press up against this switch activating the heat plate when unit is ON. It is used to automatically shut off the heat plate when the chamber is removed
- Control Knob:** This control knob adjusts the temperature of the heat plate. There are 6 levels (1-6) that can adjust the temperature from 30°C to 60°C.
- Power On Indicator (Green LED):** This LED, when illuminated, indicates that the humidifier is turned on and is ready for operation.
- Heater-On Indicator (Yellow LED):** This LED, when illuminated, indicates that the heat plate is on.
- Inlet Port with Silicon Ring:** Connects to the outlet port of the CPAP device. The Silicon Ring fits around the inlet port and ensures a tight seal between the two devices.
- Outlet port:** Connects to the flexible tubing that leads to the face/nasal mask
- AC Power Cord:** Connects to an AC outlet providing main power to the heated humidifier.
- Power Socket (Inlet):** Connects to the AC power cord from the humidifier to an outlet.
- Power Jumper Cord:** Supplies power to the CPAP by connecting to the humidifier.
- Power Socket (Outlet):** Connects to the power jumper cord
- Power Switch:** Turns the XT heated humidifier ON and OFF.

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


3. Installation


3.1 Set Up

1. Position the CPAP device onto the humidifier platform so that the four rubber knobs on the bottom of the CPAP device align with the four matching indentations on the platform. Also make sure that the Locking Tab fits squarely into the slot in the rear of the CPAP unit (beneath the power socket). This will keep the CPAP device in place on top of the humidifier.

2. Place the CPAP with humidifier onto a flat surface. This surface should be lower than the level at which you will be sleeping.

3. Connect the power jumper cord to the power socket (outlet) on the back of the humidifier and the AC inlet port on the back of the CPAP device.

⚠ NOTE: Only a 9S-005 series CPAP can be installed with this humidifier.


4. Connect the AC power cord to the power socket (inlet) on the back of the humidifier, and plug it into a main outlet.


⚠ NOTE: 1. Make sure the heated humidifier unit is suitable for the local power voltage.
2. The plug also can serve to disconnect the device.

4. Operation

4.1 Operating instruction


Daily Use

1. Turn on the humidifier by pressing the On/Off switch on the back of the humidifier platform. The humidifier power-on indicator will illuminate (green light)


2. Fill the chamber to the fill line  with distilled water (approx. 450 ml).

Silicon Ring

3. Place the chamber on the edge of the heater plate. Press down on the chamber (the heater plate will give). Slide the edges on both sides of the base of the chamber into guides on both sides of the heater plate bay. Continue to slide the chamber into place until the two air ports from the humidifier and CPAP meet and the chamber is secure against the heat plate bay of the humidifier.

Heater Plate Switch

4. Note: Make sure the water chamber is pressed up against the heater plate switch. At this time, the heater-on indicator will be lit to indicate the heater plate is on.

Note: The heater will be automatically turned off when the water chamber is removed from the heater plate.

5. Connect the flexible tubing to the outlet connector on the water chamber.

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
5. Connect the other end of the flexible tubing to your mask.

6. Adjust the heater temperature to the desired setting by using the control knob.
1 is minimum humidity, 6 is maximum humidity. The optimum setting is dependent on the ambient temperature and humidity of your room.

The Heated Humidifier and CPAP device is now ready for use. Consult the CPAP user's manual for the CPAP operating instructions.

Danger: In order to avoid damages to the CPAP equipment and prevent hazards from electric shock, DO NOT pour water into the outlet port of the water chamber while the water chamber is installed onto the humidifier. Also remember not to exceed the high water mark (Fill Line) with water when filling the chamber with water.

5. Cleaning

5.1 Cleaning the Heater

- Clean the humidifier by wiping it down with a damp cloth. Allow the heater to air dry before connecting the power cords.
- Inspect the humidifier and power cord for any damage and replace them if necessary.

5.1 Cleaning the Water Chamber

1. Turn the humidifier off and allow the heater plate and water to cool.
2. Disconnect the flexible tubing from the water chamber. Press down on the water chamber and slide it off of the heater plate. Empty the remaining water.
3. Use the extraction tool to disengage the chamber base. The extraction tool will be tucked into a recessed indentation on top of the humidifier platform. Match the notch on the front right corner of the Plastic Chamber with the inverted notch on the extraction tool so that the side of the tool is flush against the side of the chamber base. Lever the extraction tool to pry the base gently from the chamber.
4. Use a mild detergent to wash the chamber and base. Rinse all the parts with clean water and allow each part to air dry individually.
5. All items of the chamber are subject to normal wear and tear and may eventually need replacement. Replace the chamber parts if any damage is present.

6. Service

If the heated humidifier malfunctions or does not power-on, contact your equipment provider service agent immediately. Never try to open the enclosure of the heated humidifier. Maintenance for the heated humidifier requires qualified service technicians.

NOTE: In case you need to dispose of the device for any reason, you must comply with regional or national regulations for disposing such devices.

7. Troubleshooting

The table below lists troubleshooting solutions for the problems that may happen. If problems persists, contact your equipment provider service agent immediately.

Problem	Possible Causes	Solutions
Power-On indicator does not illuminate	<ol style="list-style-type: none"> 1. Power cord not plugged into a working AC outlet. 2. Power switch is not turned on 3. PCB or LED indicators are worn out 4. Fuse blown out 	<ol style="list-style-type: none"> 1. Plug the power cord into a working AC outlet. 2. See Setup procedures and turn the power on. 3. Contact your local agent or EU representative for service. 4. Contact your local agent or EU representative for service.
Heater-On indicator does not illuminate	<ol style="list-style-type: none"> 1. Heater plate is over heated 2. PCB indicator is worn out 	<ol style="list-style-type: none"> 1. Contact your local agent or EU representative for service. 2. Contact your local agent or EU representative for service.
There is no air flow through the mask	<ol style="list-style-type: none"> 1. CPAP device is not turned on or working correctly 2. Flexible tubing is not connected correctly. 3. Flexible tubing is obstructed. 	<ol style="list-style-type: none"> 1. Refer to your CPAP user manual. 2. Reconnect the flexible tubing. 3. Remove obstruction from the flexible tubing.
Condensation in mask or flexible tubing.	<ol style="list-style-type: none"> 1. The heater plate setting is too high. 2. The ambient room temperature is too cold, causing the warm humidified air in the mask to condense into droplets. 	<ol style="list-style-type: none"> 1. Adjust the control knob to lower the temperature setting. 1. Remove any air conditioner which may be in the vicinity of the heated humidifier. Or keep room temperature close to 25°C.
Water Leakage	<ol style="list-style-type: none"> 1. Assembly of water chamber is incorrect. 2. Water chamber or removable base plate worn out. 	<ol style="list-style-type: none"> 1. Remove the water chamber from the heater and discard all water. Remove base plate from the Water Chamber. Make sure the removable base plate is secured in to the base of the water chamber. Fill chamber with water until it reaches the fill line and recheck for leakage. 2. Replace the water chamber.

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
8. Technical Description

Item	Specification
Model No.:	9S-006000
Power Requirement:	120-240 VAC, 50/60Hz, 0.6-0.9A
Fuse Rating:	T2A, 250 VAC
Dimensions:	Heater: 25.5 x 17.5 x 5 cm (10" x 6.9" x 1.97") Water Chamber: 13.6 x 15 x 10.5 cm (5.36" x 5.91" x 4.14")
Weight (without Water):	Heater: under 600g (1.32 lbs) Water Chamber: under 250g (0.55 lbs)
Water Capacity:	450ml
Maximum Heater Plate Temperature:	Approximate 70°C (158°F)
Pressure Drop:	0.2cm H ₂ O @ 60LPM
Humidity Output:	10-40 mgH ₂ O/L (ambient temperature of 23°C ± 2 °C, ambient relative of humidity 60% ± 15%)
Operating Environment	Temperature: 5°C-35°C Humidity: 15%-95%RH Atmospheric Pressure: 609-768mmHg
Storage Environment	Temperature: -15°C-50°C Humidity: 10%-90%RH Atmospheric Pressure: 609-768mmHg
Classification	Class I, Type BF, IPX0 Applied Parts: Nasal Mask Not suitable for use in the presence of a flammable anesthetic mixture (No AP/APG Protection) Continuous operation.
Compliance Standard:	EN 60601-1 EN 60601-1-2 UL 60601-1 C22.2 NO. 601.1

NOTE:

1. These specifications are also suitable for other regions operating with same power supply.
2. Consult the distributor or EU representative for other technical documents.

GROUNDING INSTRUCTIONS

This product should be grounded. In the event of an electrical short circuit, grounding reduces the risk of electric shock by providing an escape wire for the electric current. This product comes with a cord equipped with a grounding wire and a grounding plug. The plug must be plugged into an outlet that is properly installed and grounded.


DANGER-Inappropriate and improper use of the grounding plug can result in a risk of electric shock.

If repair or replacement of the cord or plug is necessary, do not connect the grounding wire to either of the flat blade terminals. The insulated wire that is green (with or without yellow stripes) is the grounding wire.

Check with a qualified electrician or serviceman if you do not completely understand the grounding instructions or you are unsure if the product or outlet is properly grounded.

If it is necessary to use an extension cord, use only a 3-wire extension cord that has a three-blade grounding plug and a 3-slot receptacle that will accept the plug on the product. You can also use an adapter that has a grounding screw tab which needs to be attached to a grounded wall socket with a corresponding grounding screw. Always replace or repair a damaged cord.

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



Appendix A: EMC Information

Guidance and Manufacturer's Declaration- Electromagnetic Emissions:
This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
Harmonic emissions IEC61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

Guidance and Manufacturer's Declaration- Electromagnetic Immunity:
This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge(ESD) IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/out line	±2kV for power supply line ±1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	±1kV for differential mode ±2kV for common mode	±1kV for differential mode ±2kV for common mode	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT(>95 % dip in UT)for 0,5 cycle 40 % UT(60 % dip in UT)for 5 cycles 70 % UT(30 % dip in UT)for 25 cycles <5 % UT(>95 % dip in UT)for 5 sec	<5 % UT(>95 % dip in UT) for 0,5 cycle 40 % UT(60 % dip in UT) for 5 cycles 70 % UT(30 % dip in UT) for 25 cycles <5 % UT(>95 % dip in UT) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to the application of the test level


Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bands ^a	3Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2\sqrt{P}$ 150kHz z to 80MHz</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	<p>$d = 1.2\sqrt{P}$ 150kHz to 80MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5G MHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^c, should be less than the compliance level in each frequency range^d.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
 a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz;13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.
 b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.
 c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with

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accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device. d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0.01 0.12		0.12 0.23	
0.1 0.38		0.38	0.73
1 1.2		1.2	2.3
10 3.8		3.8	7.3
100 12		12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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