

Using the Battery with Your CPAP Machine

When using the CPAP Battery with your CPAP device, use **ONLY** the appropriate DC cord that came with your battery.

For connecting the battery pack to your CPAP machine, see specific user pages.

Turn on your machine and use as normal.

Note: Disconnecting your humidifier will increase the run time of the battery.



3

Respironics 50 Series User

Also for Puritan Bennett and DeVilbiss users.



Connect the male end of the DC cord (A) into the female end of the CPAP battery (B).

Connect the other end of the DC cord (C) into the back of your Respironics machine (D).

Respironics 60 Series User



Connect the male end of the DC cord (A) into the female end of the CPAP battery (B).

Connect the other end of your DC cord (C) into the back of your Respironics machine (D).

Please Note: The diameter of the DC Cord (C) is wider for this series.

ResMed S8 User



Connect the male end of the DC cord (A) into the female end of the CPAP battery (B).

Connect the other end of the DC cord (C) into the back of your ResMed machine (D).

6

MEDILI⁴ TECHNOLOGY LIMITED

Battery Safety and Handling continued

Do not incinerate, mutilate or expose the battery to extreme temperatures as it may cause it to burst or release toxic material.

Avoid the battery being used in an airtight compartment. The battery should be ventilated.

If the battery emits heat or shows signs of leakage, stop use.

Do not remove the outer sleeve from the battery or cut into its housing.

Disconnect the battery from your device when not in use.

During long term storage the battery should be charged and discharged once every 3 months.

Do not mix the battery with other battery brands.

Never put the battery into freshwater or seawater.

Do not attempt to take the battery apart or puncture it.

Please see page 12 for battery disposal instructions.

Battery Specifications

BATTERY PACK

Capacity: 10.4Ah

Temperature Range: 0-30°C or 32-86°F

Weight: 1.35 lbs.

CHARGER

Input: 100-240V ~ 50-60Hz 1.0A

Output: 12.6V—2.0A

Temperature Range: -10-40°C

DISPOSING OF THE BATTERY

Batteries should be discharged completely before disposal.

Do NOT dispose of the battery in your regular trash receptacle. Dispose of the battery at a facility that recycles lithium ion batteries.

Medili Technology Limited All Rights Reserved

www.medilitek.com