

Mirage Quattro™

FULL FACE MASK

FOURTH GENERATION TECHNOLOGY

Wide fit range: fourth generation full face mask from ResMed fits more than 95% of users.*

Mirage™ dual-wall cushion with spring air technology: unique cushion accommodates movement during sleep and provides an effective solution for jaw drop.

MicroFit dial: new forehead support adjustment offers personalized fitting and comfort.

Light, quiet and easy to use: these benefits are optimized with streamlined mask design, improved vent design and a range of user-friendly features.

MicroFit dial
provides better fit, seal and comfort for over 95% of users*

Set-and-forget headgear slots
save optimal fitting settings

Improved built-in vents
provide minimal noise output

Moving cushion membrane at nasal bridge
accommodates different types of nasal structures and prevents leak into the eyes

Cushion clip with tabs
allow easy removal of cushion from frame

Mirage dual-wall cushion with spring air technology
accommodates movement during sleep and provides an effective solution for jaw drop

Moving cushion membrane at chin
maintains seal and accommodates jaw drop and movement

Set-and-forget headgear clips
allow convenient removal of mask



The membrane as seen on a side view of the mask cushion

Simplified elbow with anti-asphyxia valve (AAV)
few elbow parts for easy assembly and disassembly; AAV enhances safety if the airflow to the mask is impeded





TECHNOLOGY HIGHLIGHT
MicroFit dial
 A choice of 24 positions at your fingertips—providing a comfortable, effective fit for over 95% of users*

- **Adjustable:** can be adjusted while on the user's face
- **Precise:** each turn of the dial adjusts the forehead support in small increments
- **Progressive:** moves in towards the face, allowing the user to sense when the fit is most comfortable for them
- **Works in combination with the cushion:** as the dial adjusts the angle of the mask to the face, the cushion adapts to enhance the seal where it's needed most—without the need to overtighten the headgear

USER BENEFITS

Personalized comfort

- MicroFit dial adjusts the mask to suit user's profile
- Dual-wall cushion with spring air technology adapts to facial contours and accommodates movement during sleep (including jaw drop)

Light and quiet

- Streamlined design feels light on the face
- Improved vent design reduces noise

Making it easier

- MicroFit dial and spring air cushion technology help achieve good mask fit—first time, every time
- Fewer parts for easy handling, cleaning, assembly and disassembly

PROFESSIONAL BENEFITS

Fit range increased to over 95% of users*

- MicroFit dial and spring air cushion technology adjust the mask to suit more facial structures

New XS mask size

First time fit

- First-time fit reduces callbacks and mask refitting
- Improved cushion technology seals quickly
- Simplified parts for quick and easy assembly and disassembly

Reduced inventory

- 33% fewer product codes than Ultra Mirage™ Full Face Mask
- Fewer mask parts

Comprehensive user support package

- Travel and storage bag
- Mask care DVD
- Quick-fitting guide
- User guide



PRODUCT CODES

Mirage Quattro Full Face Mask

X-Small	6120
Small	0
Medium	6120
Large	1

6120

REIMBURSEMENT CODES

2

Mirage Quattro Full Face Mask 6120

3

HCPCS Code Descriptor

A7030	Fullface mask used with positive airway pressure (PAP) device, 1 per 3 months
A7031	Replacement cushion for full face mask, 1 per month
A7035	Headgear used with positive airway pressure (PAP) device, 1 per 6 months

* Source: internal ResMed anthropometric data