

AirSense™ 11 CPAP devices

At-home or hospital PAP therapy that treats obstructive sleep apnea (OSA) in patients weighing more than 66 lbs (30 kg).



	AirSense 11 AutoSet™ Advanced self-adjusting AutoSet therapy device	AirSense 11 Elite Premium CPAP therapy device	AirSense 11 CPAP Basic CPAP therapy device
Operating modes and pressures			
CPAP	•	•	•
AutoSet or AutoSet for Her	•	–	–
Operating pressure range	4–20 cm H ₂ O	4–20 cm H ₂ O	4–20 cm H ₂ O
AutoSet mode default pressure	Min: 4 / Max: 20	–	–
Features			
AutoSet Response comfort setting	•	–	–
CSA, RERA and CSR detection	•	•	–
AutoRamp with sleep onset detection	•	•	•
Easy-Breathe™ waveform with EPR	•	•	•
SmartStart™/Stop feature	Individually toggleable	Individually toggleable	Individually toggleable
Oxygen use	Up to 15 L/min	Up to 15 L/min	Up to 15 L/min
Remote diagnostic information in AirView	•	•	•
Patient access to ResMed myAir™	•	•	•
Bluetooth® connectivity	•	•	•
Over-the-air updates	•	•	•
Humidification			
Standard heated humidifier	Integrated	Integrated	Integrated
ClimateLineAir™ 11 heated tubing	Optional	Optional	Optional
Data transfer on SD card			
Detailed data	30 sessions	30 sessions	–
High-resolution flow	30 sessions	30 sessions	–
Compliance and summary data	365 sessions	365 sessions	365 sessions

Model	Part number (US)	Part number (Canada)	HCPC	Recommended replacement schedule [†]
AirSense 11 AutoSet				
• With HumidAir™	39000	39006	E0601	Every 5 years
• With HumidAir and ClimateLineAir	39001	39007		
AirSense 11 Elite				
• With HumidAir	39004	39010	E0601	Every 5 years
• With HumidAir and ClimateLineAir	39005	39011		
AirSense 11 CPAP				
• With HumidAir	39002	39008	E0601	Every 5 years
• With HumidAir and ClimateLineAir	39003	39009		

[†]The myAir by ResMed app is available in English and Spanish in the US and in English and French in Canada.

[†]Replacement cycle may vary depending on patient's insurance coverage.

Air11 accessories

	Accessory	Part number	HCPC	Recommended replacement schedule*
Humidification				
	HumidAir™ 11 standard tub	39100	A7046	Every 6 months
	HumidAir 11 cleanable tub	39101	A7046	Every 6 months
	ClimateLineAir™ 11 heated tubing	39102	A4604	Every 3 months
	SlimLine™ tubing	36810	A7037	Every 3 months
	6' standard tubing	14994	A7037	Every 3 months
	6'6" standard tubing	14987	-	-
	6'6" autoclavable standard tubing	14986		
3m tubing	14999			
	Air11 side cover	39226	-	-
	Air11 air outlet	39220		
Filter				
	Air11 standard filters, 1 pk	39300	A7038	2 per month
	Air11 standard filters, 2 pk	39301		
	Air11 standard filters, 12 pk	39302		
	Air11 standard filters, 50 pk	39303		
	Air11 hypoallergenic filters, 1 pk	39304		
	Air11 hypoallergenic filters, 2 pk	39305		
	Air11 hypoallergenic filters, 12 pk	39306		
	Air11 hypoallergenic filters, 50 pk	39307		
	Filter door replacement	29108	-	-
Power,  connectivity				
	AirSense™ 11 power supply unit, 65W	39205	-	-
	AirCurve™ 11 power supply unit, 90W	39206		
	Air11 DC/DC converter	39231		
	Air11 SD card with envelope, 1 pk	39228		
	AirSense 11 SD card with envelope, 10 pk	39229		
	Air11 SD card door	39225		
	AirCurve 11 SD card door	39233		
	Air11 BLE WristOx-2	7079666		
	Air11 expandable wrist band	7079667		
	Air11 small soft sensor, W02 connector	7079814		
	Air11 medium soft sensor, W02 connector	7079815		
	Air11 single-use wrist band, 15 pk	7079668		
	Air11 large soft sensor, W02 connector	7079816		
Air11 single-use sensor and wrist band, 15 pk	70710011			
Faceplates,  quick start guides				
	Air 11 faceplate	29106	-	-
	Air11 travel bag	39221		
	AirSense 11 user guide	398133		
	AirCurve 11 user guide	398195		
	Air11 quick start guide	398122		

*Replacement cycle may vary depending on patient's insurance coverage.