

AirSense™ 10 CPAP devices

For at-home or hospital treatment of obstructive sleep apnea (OSA) in patients weighing more than 66 lbs (30 kg).



	AirSense 10 AutoSet™ Advanced self-adjusting AutoSet therapy device	AirSense 10 Elite Premium CPAP therapy device	AirSense 10 CPAP Basic CPAP therapy device
Operating modes and pressures			
CPAP	•	•	•
AutoSet	•	–	–
Operating pressure range	4–20 cm H2O	4–20 cm H2O	4–20 cm H2O
AutoSet mode default pressure	Min: 4 / Max: 20	–	–
Features			
AutoSet Response comfort setting	•	–	–
CSA, RERA and CSR detection	•	•	–
AutoRamp with sleep onset detection	•	•	•
Easy-Breathe™ waveform with EPR	•	•	•
SmartStart™/Stop feature	•	•	•
Remote diagnostic information in ResMed AirView™	•	•	•
Integrated cellular module	•	•	•
Patient access to ResMed myAir™†	•	•	•
Oxygen use	Up to 4 L/min	Up to 4 L/min	Up to 4 L/min
Standard heated humidifier			
ClimateLineAir™ 10 heated tubing	Integrated	Integrated	Integrated
	Optional	Optional	Optional
Detailed data			
High-resolution flow	30 sessions	30 sessions	–
Compliance and summary data	30 sessions	30 sessions	–
	365 sessions	365 sessions	365 sessions

Model	Part number (US)	Part number (Canada)	HCP#	Recommended replacement schedule‡
AirSense 10 AutoSet				
• With HumidAir™	3720	37402	E0601	Every 5 years
• With HumidAir and ClimateLineAir	7	37403		
AirSense 10 Elite	3720			
• With HumidAir	8	3740	E0601	Every 5 years
• With HumidAir and ClimateLineAir	3720	0		
AirSense 10 CPAP	3720	3740		
• With HumidAir	6	1	E0601	Every 5 years
• With HumidAir	37203	37398		
• With HumidAir and ClimateLineAir	37204	37399		

† ASV therapy is contraindicated in patients with chronic, symptomatic heart failure (NYHA 2-4) with reduced left ventricular ejection fraction (LVEF ≤ 45%) and moderate to severe predominant central sleep apnea.

‡ The myAir by ResMed app is available in English and Spanish in the US and in English and French in Canada.

§ Replacement cycle may vary depending on patient's insurance coverage.

Air10 accessories

	Accessory	Part number	HCPC	Recommended replacement schedule*
Humidification				
	HumidAir™ standard tub HumidAir	37299	A7046	Every 6 months
	cleanable tub	37479	A7046	Every 6 months
	ClimateLineAir™ heated tubing SlimLine™ tubing 6" standard tubing 6" standard tubing 6" autoclavable standard tubing	37298 36810 14994	A4604 A7037 A7037	Every 3 months Every 3 months Every 3 months
	3m tubing HumidAir side cover, gray (AirSense™ 10 AutoSet for Her, AirCurve™ 10 ST)	14987 14986 14999		
	HumidAir side cover, charcoal (AirSense 10 AutoSet and Elite, AirCurve 10 ST-A)	37335	-	-
		37303		
Filters				
	Air10 standard filters, 1 pk	36850	A7038	2 per month
	Air10 standard filters, 2 pk	36851		
	Air10 standard filters, 12 pk	36852		
	Air10 standard filters, 50 pk	36853		
	Air10 hypoallergenic filters, 1 pk	36855		
	Air10 hypoallergenic filters, 2 pk	36856		
	Air10 hypoallergenic filters, 12 pk	36857		
	Air10 hypoallergenic filters, 50 pk	36858		
	Filter cover, gray (AirSense 10 AutoSet for Her, AirCurve 10 ST)	19768		
	Filter cover, charcoal (AirSense 10 AutoSet and Elite, AirCurve 10 ST-A)	19729	-	-
Power, data management and oximetry				
	Air 10 power supply unit	37474	-	-
	Air10 DC/DC converter	37297		
	RPS II-V2 kit for Air10	24702		
	RPS II-V2 battery	24705		
	RPS II DC cable for Air10	37343		
	RPS II-V2 adapter for Air10	24706		
	Air10 SD card protective folder, with card, 1 pk	37378		
	Side panel SD card door, gray (AirSense 10 AutoSet for Her, AirCurve 10 ST)	19766		
	Side panel SD card door, charcoal (AirSense 10 AutoSet and Elite, AirCurve 10 ST-A)	19728		
	Air10 oximetry kit: Module; XPOD; XPOD clip and reusable soft sensor (medium)	370004		
	Air10 oximetry kit module	37302		
	XPOD LP oximeter	22374		
	XPOD LP clip	22371		
	Reusable soft sensor (medium)	707560		
Faceplates, travel bags and quick start guides				
	Air10 faceplate, silver	19736	-	-
	Air10 faceplate, white with gray	19735		
	Air10 faceplate, black	19734		
	Air10 travel bag	37304		
	Air10 travel bag for her	37305		

*Replacement cycle may vary depending on patient's insurance coverage.