



## Mask Frame Sizing


- Medium (MED) frame - 20" in length
- Small (SM) frame - 19" in length
- Large (LG) frame - 22" in length



The medium (MED) mask frame will comfortably fit most faces. If the MED frame does not fit your face, contact your healthcare provider to see if a small (SM) or large (LG) mask frame would better suit your needs.



\* **Tip:** If the mask frame falls backward and is too close to your ears you may need a smaller mask frame.



\* **Tip:** If the mask frame falls forward on the head and is too close to your eyes you may need a larger mask frame.