



**ResMed**

Actual Size: 5 1/2" (Width)  
x 6 7/8" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

**AirFit™ F30**

Full face mask



© 2018 ResMed Ltd. All Rights Reserved.  
Specifications may change without notice.  
Product availability may vary across regions.  
10111032/1 2018-10

**Rx Only** CE 0123

