

Mirage Micro™

NASAL MASK

MICRO MAKES (F)IT EASY

The latest generation nasal mask from ResMed provides personalized fit, seal and comfort for more users.

Combined with ResMed's proven Mirage™ cushion technology, the Mirage Micro's unique MicroFit dial fine-tunes fit to suit each individual user. Intuitive to use for sleep professionals and users,

the Mirage Micro includes two cushion sizes packaged together (except small) for fewer product codes and reduced inventory. With enhanced comfort and an improved fit range of over 95%*, the Mirage Micro is the new benchmark in nasal masks.

MicroFit dial

Offers 24 positions for better fit, seal and comfort

Streamlined forehead support

Provides clear field of vision due to sleek, stylish design

Quiet vent design

Reduces noise output; includes removable clip for easy cleaning

360° elbow rotation

Multiple tubing positions for maximum comfort

Set-and-forget headgear clips

Retains optimal headgear settings; easy to attach and detach

Elbow partition

Aids CO₂ washout and reduces noise output

Cushion loop

For easy removal of cushion


Mirage dual-wall cushion

Reduces pressure on nasal bridge and seals effectively; clipless for easy assembly and cleaning



Swivel

Offers a convenient way to connect and disconnect air tubing





TECHNOLOGY HIGHLIGHT

MicroFit dial

A choice of 24 positions at your fingertips—providing a comfortable, effective fit for more users

- **Adjustable:** can easily be adjusted while on the user's face
- **Precise:** each turn of the dial adjusts the forehead support in small increments
- **Progressive:** moves in toward the face, allowing the user to sense when the fit is most comfortable for them
- **Works in combination with the cushion:** as the dial adjusts the angle of the mask to the face, the dual-wall cushion adapts to enhance the seal without the need for over-tightening the headgear

USER BENEFITS

Personalized fit and comfort

- MicroFit dial fine-tunes fit to suit each user's profile
- Mirage dual-wall cushion is available in five sizes, including new Small size, to fit more users
- Deep mask frame design accommodates a greater range of nasal structures

Quiet and streamlined

- Slender design for improved field of vision
- Vent design quietly disperses air

Intuitively easy

- MicroFit dial is easily adjusted with one hand, even while lying down
- Fewer mask parts for easy cleaning & assembly
- Comprehensive user support package includes travel bag and user guide with mask fitting and care instructions

PROFESSIONAL BENEFITS

First-time fit for more users

- Combination of MicroFit dial and dual-wall cushion achieves a perfect fit for a range of facial structures—first time, every time
- Five cushion sizes fit over 95% of users*

Customer satisfaction and compliance

- Increased therapy acceptance through first-time fit
- Fewer mask parts for fast and intuitive assembly and disassembly

Cost and time efficiencies

- Reduced inventory: three product codes for five cushion sizes
- Reduced fitting and training time
- Reduced callbacks: ease of use and fitting



PRODUCT CODES

Mirage Micro Nasal Mask

Americas

Small	16333
Medium & Large	16334
Large Wide & X-Large	16335

REIMBURSEMENT CODES

Mirage Micro Nasal Mask

HCPCS	Code Descriptor
A7034	Nasal mask used with positive airway pressure device, 1 per 3 months
A7032	Replacement cushion for nasal mask, 2 per month
A7035	Headgear used with positive airway pressure device, 1 per 6 months

*Based on ResMed anthropometric database