



The right fit for you



**AirFit™** F20

Full face mask

ResMed.com

# Start dreaming again

The benefits of sleep apnea treatment are truly remarkable. Getting a good night's sleep gives you the energy to get back to doing the things you love. But one of the most inspiring things people tell us about starting CPAP is that they're able to dream again.





Congratulations on taking an important step toward getting restful sleep.

ResMed is the most widely-used and #1 mask brand preferred by patients<sup>1</sup>, and we're excited to present the AirFit™ 20 series.

With feedback from more than 350 patients worldwide used to develop the AirFit 20 series<sup>2-5</sup>, you can feel confident that our most technologically advanced full face mask yet is also the right fit for you.

**AirFit™** F20

Full face mask




# Ultimate comfort and performance

The AirFit F20™ features an InfinitySeal™ cushion – the most adaptive we've ever designed – along with a soft and flexible frame that includes integrated padding. The result is a mask that comfortably moves with you and provides a reliable seal throughout the night, regardless of sleep position.

- 1 **InfinitySeal™ cushion** is proven to fit 97% of all faces.<sup>2-4</sup>
- 2 **Quick-release elbow** makes it simple to disconnect from the mask tubing without removing the mask.
- 3 **Plush headgear** is specifically designed for extra softness and stability.
- 4 **Soft and flexible frame** sits underneath the eyes and adapts to different facial types, while its integrated padding provides added softness.
- 5 **Magnetic clips** guide the headgear to the frame quickly and easily.







# Sleep support tailored to you

- Register your mask
- Get reminders about part replacement
- Receive personalized treatment tips

Get started at  
[ResMed.com/RegisterF20](https://ResMed.com/RegisterF20)



**1** Combined responses of those expressing a preference in an independent patient survey in the UK, USA, France, Germany and Australia between April 10 to May 11, 2015. Visit [ResMed.com/MaskBrand](https://ResMed.com/MaskBrand). **2** ResMed AirFit F20 internal Australian fitting study of existing CPAP patients conducted between March–April 2016. **3** ResMed AirFit F20 internal USA fitting study of existing CPAP patients conducted April 2016. **4** ResMed AirFit F20 internal EU and APAC fitting study of existing CPAP patients conducted June 2016. **5** ResMed AirFit N20 internal international fitting study of new and existing CPAP patients, conducted Nov 2015.

Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA  
+1 858 836 5000 or 1 800 424 0737 (toll free). See [ResMed.com](https://ResMed.com) for other ResMed locations worldwide. AirFit and InfinitySeal are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. © 2016 ResMed. 1019457/1 2016-10

[ResMed.com](https://ResMed.com)

Follow ResMed on:

