



## AirMini™

AirMini is a self-adjusting system indicated for the treatment of obstructive sleep apnea (OSA) in patients, including female patients with mild to moderate OSA when using AutoSet™ for Her treatment mode, weighing more than 66 lb (30 kg). Its small size makes it an ideal portable PAP, and the AirMini app allows patients to control a range of comfort settings.



Model	Part number (US)
AirMini	3 8
AirMini mask pack N30 (mask included)	113
AirMini mask pack P10 (mask included)	38878
AirMini setup pack N20 (mask not included)	38824
AirMini setup pack F20 (mask not included)	38823
AirMini setup pack F30 (mask not included)	38011

	Accessory	Part number
<b>For nasal and nasal pillows masks</b>		
	HumidX, 3 pk	38809
	HumidX, 6 pk	38810
	HumidX, 50 pk	38808
	HumidX Plus, 3 pk	38812
	HumidX Plus, 6 pk	38813
	HumidX Plus, 50 pk	38811
<b>For full masks</b>		
	HumidX F20, 3 pk	38012
	HumidX F20, 6 pk	38013
<b>Travel</b>		
	AirMini travel bag	38840
	AirMini travel case	38841
<b>Filters</b>		
	AirMini filter, 1 pk	38814
	AirMini filter, 2 pk	38815
	AirMini filter, 12 pk	38816
	AirMini filter, 50 pk	38837
<b>Mount system</b>		
	AirMini mount system	38842
	<b>Power</b>	
	AirMini DC/DC converter	38839
	AirMini 20W AC adapter	38828

\*Not available with the AirFit™ F30 mask.