

## DreamStation Sleep Therapy System Configurations:

	DreamStation CPAP	DreamStation CPAP Pro	DreamStation Auto CPAP	DreamStation BiPAP Pro	DreamStation Auto BiPAP
<b>Modes</b>	Fixed CPAP	Fixed CPAP CPAP-Check	Fixed CPAP CPAP-Check Auto-CPAP	Fixed CPAP Fixed BiPAP	Fixed CPAP Fixed BiPAP Auto-BiPAP
<b>Features</b>	n/a	EZ-Start Auto-Trial	EZ-Start Auto-Trial Opti-Start	EZ-Start	EZ-Start
<b>Pressure Range</b>	4 to 20 cm H <sub>2</sub> O	4 to 20 cm H <sub>2</sub> O	4 to 20 cm H <sub>2</sub> O	4 to 25 cm H <sub>2</sub> O	4 to 25 cm H <sub>2</sub> O
<b>Flex Comfort</b>	C-Flex	C-Flex C-Flex+	C-Flex C-Flex+ A-Flex	C-Flex C-Flex+ Bi-Flex	C-Flex C-Flex+ Bi-Flex
<b>Humidification Modes (3)</b>	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube	Fixed Adaptive Heated Tube
<b>Mask Type Resistance Control</b>	Standard	Standard	Standard	Standard	Standard
<b>Ramp (0-45 minutes)</b>	Ramp SmartRamp	Ramp SmartRamp	Ramp SmartRamp	Ramp SmartRamp	Ramp SmartRamp
<b>Altitude Compensation</b>	Automatic	Automatic	Automatic	Automatic	Automatic
<b>Advanced Event Detection</b>	n/a	Standard	Standard	Standard	Standard
<b>Flow Waveforms</b>	n/a	Standard	Standard	Standard	Standard
<b>Modem Compatibility</b>	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)	Standard (optional Cellular Modem or Wi-Fi Accessory)
<b>Built-in Bluetooth Connectivity</b>	Standard	Standard	Standard	Standard	Standard
<b>DreamMapper Compatibility</b>	Standard (by mobile app or web)	Standard (by mobile app or web)	Standard (by mobile app or web)	Standard (by mobile app or web)	Standard (by mobile app or web)
<b>Oximetry Compatibility</b>	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)	Standard (optional Nonin accessories)
<b>Check Mask Fit feature</b>	Standard	Standard	Standard	Standard	Standard
<b>Performance Check</b>	Standard	Standard	Standard	Standard	Standard