



7 Mask fitting

To choose a well fitted mask, ensure the entire mouth and nose are included inside the mask.

Make sure the top of the cushion is not too high on the bridge of the nose area. No leak should remain.

Landmarks

- a** Side of mouth
- b** Bridge of nose zone
- c** Below the lower lip



Mask size (inches)	Length	Width	Height
Infant	2.76	2.60	1.93
Toddler	3.46	3.15	2.32
Small adult	4.53	3.74	3.19
Medium adult	4.84	4.02	3.54
Large adult	4.92	4.33	3.54



Mouthpiece



Trach adaptor



Toddler mask



Small, medium and large masks