



## Exhaust-flow information

F&P Vitera™		F&P Simplus™		Forma™		FlexiFit™432		FlexiFit 431		F&P Brevida™		F&P Pilairo™ Q		Opus™ 360	
Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)**	Pressure (cmH <sub>2</sub> O)	Flow (l/min)**	Pressure (cmH <sub>2</sub> O)	Flow (l/min)
4	24	4	22.5	4	20.5	4	22.5	4	26	4	21	4	18	4	19.5
5	26	5	25	5	23	5	24	5	28	5	24	5	20	5	22
6	28	6	27	6	24.5	6	26	6	30	6	27	6	22	6	24
7	30	7	29.5	7	26	7	28	7	32	7	29	7	24	7	26
8	32	8	31.5	8	27.5	8	29.5	8	33.5	8	31	8	26	8	28
9	34	9	33.5	9	29	9	31	9	35	9	33	9	28	9	30
10	36	10	35	10	29.5	10	33.5	10	36	10	35	10	30	10	31.5
11	38	11	37	11	30	11	36	11	37	11	37	11	31	11	33
12	40	12	38.5	12	31	12	37.5	12	38.5	12	39	12	33	12	34.5
13	42	13	40.5	13	32	13	39	13	40	13	41	13	34	13	36
14	43	14	42	14	33.5	14	40	14	41.5	14	43	14	36	14	37.5
15	45	15	43.5	15	35	15	41	15	43	15	44	15	38	15	39
16	47	16	45	16	36	16	42.5	16	43.5	16	46	16	39	16	40.5
17	48	17	46.5	17	37	17	44	17	44	17	47	17	41	17	42
18	50	18	47.5	18	38.5	18	45.5	18	45	18	49	18	42	18	43.5
19	51	19	49	19	40	19	47	19	46	19	50	19	44	19	45
20*	53	20	50.5	20	42.5	20	48.5	20	47	20	52	20	45	20	46.5

  

F&P Evora™		F&P Eson™ 2		F&P Eson™		F&P Zest™Q		Zest™ Range		FlexiFit™407		FlexiFit 406		FlexiFit 405	
Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)**	Pressure (cmH <sub>2</sub> O)	Flow (l/min)**	Pressure (cmH <sub>2</sub> O)	Flow (l/min)**	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)	Pressure (cmH <sub>2</sub> O)	Flow (l/min)**
4	17	4	19	4	20	4	18	4	20	4	22	4	21	4	21
5	19	5	22	5	22	5	20	5	23	5	25	5	24	5	24
6	21	6	25	6	25	6	22	6	25.5	6	27.5	6	26	6	26.5
7	23	7	27	7	27	7	24	7	28	7	30	7	28	7	29
8	24	8	29	8	29	8	26	8	30	8	32	8	30	8	31
9	26	9	32	9	31	9	28	9	32	9	34	9	32	9	33
10	28	10	34	10	33	10	29	10	34	10	36	10	34	10	35
11	29	11	36	11	34	11	31	11	36	11	38	11	36	11	37
12	31	12	38	12	36	12	32	12	38	12	40	12	37.5	12	39
13	32	13	40	13	38	13	33	13	40	13	42	13	39	13	41
14	34	14	41	14	39	14	35	14	42	14	43.5	14	40.5	14	42.5
15	35	15	43	15	41	15	37	15	44	15	45	15	42	15	44
16	37	16	45	16	43	16	38	16	45.5	16	46.5	16	43.5	16	45.5
17	38	17	47	17	44	17	39	17	47	17	48	17	45	17	47
18	39	18	48	18	45	18	40	18	48	18	49.5	18	46.5	18	48.5
19	40	19	50	19	47	19	41	19	49	19	51	19	48	19	50
20***	41	20	52	20	48	20	42	20	50.5	20	52.5	20	49	20	51.5

\*F&P Vitera operates up to 30 cm H<sub>2</sub>O. To view the complete exhaust flow table refer to the F&P Vitera Use & Care Guide.  
 \*\*When used with Diffuser and Cover.  
 \*\*\*F&P Evora operates up to 25 cm H<sub>2</sub>O. To view the complete exhaust flow table refer to the F&P Evora Use & Care Guide.