

## Aloha™ Nasal Pillow System

### Lightweight Comfort

The Aloha Nasal Pillows System is the next evolution in the design of the nasal pillow interface. Not only is this device **extremely lightweight and quiet**, but it offers an **incredibly stable and comfortable seal**.

The unique ball-and-swivel joint, along with a hose swivel, **adapts to the most active of sleepers** – all with minimal contact to the user's face.

The Aloha features a unique and exclusive system that allows the end-user to adjust the angle and depth of the nasal pillows – improving comfort and seal while minimizing potential discomfort.



Aloha Nasal  
Pillow System  
(ALO125)

### Performance Features

- Ball-and-swivel joint minimizes torque and pull to accommodate various sleeping positions
- **Exclusive Arced-Track™ technology** allows the angle of the pillow reservoir to be adjusted for a more customized and natural fit
- CPAP hose can be routed over the head or worn down the chest
- Wide-set headgear provides limited points of contact on the face and no interference with the eyes
- Angled exhalation ports minimize potential disruption or discomfort for the user or bed partner
- Minimalistic headgear made of soft, pliable, and breathable material



## Design Features

- Latex-Free
- Four-Point adjustable headgear
- Removable and ultra-soft micro-fleece headgear wraps
- Simple to assemble and disassemble – one piece reservoir design
- Supple, soft and stable nasal inserts
- Tubing Ring to secure hose, when routed over-the-head to headgear
- Headgear wraps included to maximize comfort

## Reimbursement Information

### Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
Aloha Pillows	Pillows for nasal application device	A7033	Two pairs per month	\$22.36 - \$26.30
Aloha Application	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	One per 3 months	\$92.57 - \$108.90
Headgear	Headgear used with positive airway pressure device	A7035	One per 6 months	\$31.27 - \$36.79

The noted allowable rates are based off of the January 2012 HCPCS Fee Schedule. All information presented is for reference purposes only.

## Ordering Information

ITEM NUMBER	PRODUCT DESCRIPTION	ITEM NUMBER	PRODUCT DESCRIPTION
AL0100	Complete Aloha Nasal Pillow System (All Size Kit)	AL0130	Aloha Headgear
AL0121	Aloha Pillows – Small	AL0131	Aloha Frame – Small
AL0123	Aloha Pillows – Medium	AL0133	Aloha Frame – Medium
AL0125	Aloha Pillows – Large	AL0135	Aloha Frame – Large



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
 800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)