

Aloha™ Nasal Pillow System

Lightweight Comfort

The Aloha Nasal Pillows System is the next evolution in the design of the nasal pillow interface. Not only is this device **extremely lightweight and quiet**, but it offers an **incredibly stable and comfortable seal**.

The unique ball-and-swivel joint, along with a hose swivel, **adapts to the most active of sleepers** – all with minimal contact to the user's face.

The Aloha features a unique and exclusive system that allows the end-user to adjust the angle and depth of the nasal pillows – improving comfort and seal while minimizing potential discomfort.



Aloha Nasal
Pillow System
(ALO125)

Performance Features

- Ball-and-swivel joint minimizes torque and pull to accommodate various sleeping positions
- **Exclusive Arced-Track™ technology** allows the angle of the pillow reservoir to be adjusted for a more customized and natural fit
- CPAP hose can be routed over the head or worn down the chest
- Wide-set headgear provides limited points of contact on the face and no interference with the eyes
- Angled exhalation ports minimize potential disruption or discomfort for the user or bed partner
- Minimalistic headgear made of soft, pliable, and breathable material



Design Features

- Latex-Free
- Four-Point adjustable headgear
- Removable and ultra-soft micro-fleece headgear wraps
- Simple to assemble and disassemble – one piece reservoir design
- Supple, soft and stable nasal inserts
- Tubing Ring to secure hose, when routed over-the-head to headgear
- Headgear wraps included to maximize comfort

Reimbursement Information

Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
Aloha Pillows	Pillows for nasal application device	A7033	Two pairs per month	\$22.36 - \$26.30
Aloha Application	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	One per 3 months	\$92.57 - \$108.90
Headgear	Headgear used with positive airway pressure device	A7035	One per 6 months	\$31.27 - \$36.79

The noted allowable rates are based off of the January 2012 HCPCS Fee Schedule. All information presented is for reference purposes only.

Ordering Information

ITEM NUMBER	PRODUCT DESCRIPTION	ITEM NUMBER	PRODUCT DESCRIPTION
AL0100	Complete Aloha Nasal Pillow System (All Size Kit)	AL0130	Aloha Headgear
AL0121	Aloha Pillows – Small	AL0131	Aloha Frame – Small
AL0123	Aloha Pillows – Medium	AL0133	Aloha Frame – Medium
AL0125	Aloha Pillows – Large	AL0135	Aloha Frame – Large



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
 800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com