





DreamStation pulse oximetry information

Setup of oximetry requires 3 components:

1. DreamStation Cellular Modem, US with oximetry; DreamStation Wi-Fi Accessory, US with oximetry; or a DreamStation Link Module
2. Nonin® Xpod SpO₂ assembly
3. Nonin oximetry sensors

Once components are connected:

- DreamStation will display the “good connection” icon  if the PAP oximetry sensor is connected to the patient, working, and therapy is on
- DreamStation displays the “Good SpO₂” icon  when the PAP is turned off, after 4 or more hours of therapy and oximetry data has been collected



| Description | Part number |
|--|-------------|
| DreamStation Cellular Modem, US with oximetry | 100600C |
| DreamStation Wi-Fi Accessory, US with oximetry | 100700W |
| DreamStation Link Module | 1120293 |
| Nonin SpO ₂ assembly | 1121694 |
| Reusable finger sensor (clip) | 936 |
| Oximetry finger sensor, adult flex | 953 |
| Flexiwrap, sensor tape, adult, 25/pk | 954A |
| Reusable finger sensor (soft), medium | 1062050 |
| Reusable finger sensor (soft), small | 1062051 |

Note: DreamStation is only compatible with Nonin pulse oximetry components