



**CPAP Battery Pack Information:**

Power :115W

**The CPAP Battery Pack comes complete with:**

CPAP Battery + Power Cord + Charger + User Guide+ Carry Bag

**Estimated Run Times (In Hours) for Various Machine Models:**

CPAP(Average Usage Time in Hours)

Pressure (cm H2O)	Respironics REMstar / REMstar M Series	Philips Respironics System One REMstar	ResMed S8	ResMed S9
6	20	24	17	17
7	18	19	15	15
8	17	16	14	13
9	16	14	13	12
10	14	12	12	11
11	13	11	11	11
12	11	11	10	10
13	10	10	9	9
14	9	10	9	9
15	8	9	8	8
16+ **	<8	<9	<8	<8

Note that these run times are approximate and depend on several factors, including CPAP machine type and inverter type (if an inverter is being used). The use of a heater on a humidifier will substantially reduce battery run time, probably by more than 50%.