

# Humidifier/device compatibility list

English

This supplement specifies the devices compatible with ResMed's humidifiers. For full details regarding the correct use of these devices, please refer to the relevant User Guides.

Note: Not all products are available in all regions.

Device	Heated Humidifiers						External Humidifier	Passive Humidifiers	
	HumidAire 2i™	HumidAire 3i™ H3i Plus™	H4i™	H4i™ for Stellar	H5i™ / H5i™ for Her	HumidAir™	HumiCare™ D900	Passover	HumidAire 2iC™
AirSense™ 10						✓			
AutoSet Spirit™	✓							✓	✓
S7™ Series	✓							✓	✓
S8 AutoSet Vantage™		✓	✓						
S8 Escape™		✓	✓					✓	
S8 Escape™ II		✓	✓						
S8 Elite™ II, S8 AutoSet™ II		✓	✓						
S9 AutoSet™, S9 AutoSet™ for Her, S9 Elite™, S9 Escape™, S9 Escape™ Auto					✓				
VPAP™ III Series	✓							✓	✓
VPAP™ III ST-A with QuickNav	✓								✓
VPAP™ Adapt SV	✓								✓
VPAP™ Adapt (S9)					✓				
VPAP™ Adapt (S9 ASV)					✓		✓		
VPAP™ Auto		✓	✓						
VPAP™ Auto (S9)					✓				
VPAP™ S (S9), VPAP™ ST (S9), VPAP™ ST-A (S9 iVAPS), VPAP™ COPD (S9)					✓		✓		
VPAP™ Auto 25, VPAP™ ST			✓						
VPAP™ Tx	✓								
S9 VPAP™ Tx					✓				
VPAP Malibu™	✓								
Lumis™ Tx					✓		✓		
Stellar™				✓			✓		



Stellar **ResMed Germany Inc.** Fraunhoferstr. 16  
82152 Martinsried Germany

All others **ResMed Ltd** 1 Elizabeth Macarthur Drive Bella Vista  
NSW 2153 Australia

#### DISTRIBUTED BY

**ResMed Corp** 9001 Spectrum Center Boulevard San Diego CA 92123 USA  
**ResMed Ltd** 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

**ResMed Corp** 9001 Spectrum Center Boulevard San Diego CA 92123 USA  
**[EC/REP] ResMed (UK)** (for Ltd) 96 Jubilee Ave Milton Park Abingdon Oxfordshire OX14 4RW UK

See [www.resmed.com](http://www.resmed.com) for other ResMed locations worldwide. For patent information, see [www.resmed.com/ip](http://www.resmed.com/ip). AirSense, HumidAir, S8, S9, Lumis Tx, Stellar, PaceWave and VPAP are trademarks of ResMed Ltd. S9, Stellar, and VPAP are registered in U.S. Patent and Trademark Office. © 2014 ResMed Ltd. 1017966/1 2014-08.