

## SoClean 2



ideal for at-home use

## SoClean 2 Go



ideal for travel

Fully Sanitizes CPAP Mask, Hose & Reservoir

YES

YES

Works with a Wide Variety of Home CPAP Machines

YES

YES

Works with a Wide Variety of CPAP Masks

YES

YES

## SoClean 2



ideal for at-home use

## SoClean 2 Go



ideal for travel

Works with Heated Hoses (with Adapter)

YES

YES

Ideal for Travel

NO

YES

Programmable Timer

YES

NO

Status Indicator Lights

YES

YES

Automatic Shut-Off

YES

YES

CPAP Hose Entry from Left or Right

YES

N/A

## SoClean 2



ideal for at-home use

## SoClean 2 Go



ideal for travel

Check Valve Assembly

Replace Every Six Months

Replace Every Six Months

Filter

Replace Every Six Months

Replace Every Six Months

Outer Dimensions

7.5"L x 6.5"W x 8.5"H

4.3"L x 3.9"W x 4.9"H

Footprint Dimensions

7.5"L x 6.5"W

4.3"L x 3.9"W

Chamber Dimensions

5.6"L x 5.4"W x 5"D

N/A (Bag)

Chamber Volume

178.3 in<sup>3</sup>

N/A

## SoClean 2



ideal for at-home use

## SoClean 2 Go



ideal for travel

Bag Dimensions

N/A

10.7"L x 10.5"W

Bag Volume

N/A

250 in<sup>3</sup>

Weight

5.6 lbs

1.3 lbs

Cord Length

4' 8"

cordless

Electrical

12 VDC

12 VDC

Power Supply

110-240 VAC

Battery Operated

## SoClean 2



ideal for at-home use

## SoClean 2 Go



ideal for travel

Power Consumption

5 watts

5 watts

Warranty

1 year

1 year