



# ZOPEC EXPLORE Oxygen

*Professional Sleep Guardian*



**576 Wh with Platinum UPS Technology**

## Zopec Explore Oxygen Universal CPAP, BiPap, and Oxygen Battery (the Beast!)

- **Common Uses:** UPS protection using CPAP or BiPap with humidifier and heated tube.
- **Platinum Level UPS Technology:** Uninterrupted sleep when power outages occur!
  - \* This UPS supports running humidifier and heated tube (i.e. ResMed AutoClimate Control)
  - \* Battery stays at 100% as standby power. Auto-Switch to battery power instantly with sensors.
  - \* Doesn't hurt battery life plugged in 24/7! Normally uses wall power with pass-through.
- **Mega Battery Capacity:** 576 Wh / 160 Ah. (about 2 car batteries but only 7.5 lb.!).
  - \* Sleep 8 - 12 hours with humidifier and heated tube. (i.e. ResMed AutoClimate Control)
  - \* Sleep 16 - 24 hours with humidifier (medium setting) and no heated tube
  - \* Sleep 42 - 64 hours without humidifier and heated tube (up to 8 nights of sleep!)
- **Solid State Peltier Cooling:** Battery is 100% silent! No cooling fans.
- **Panasonic Lithium Ion Cells:** Same brand as Tesla vehicles use! Smallest, lightest, best quality!
- **Super Light and Compact:** 7.5 lb. 10.25" x 6.25" x 3.25" in. (2 car batteries weigh 80 lb.!).
- **Two 110V AC Outlets (world):** Works with all machine brands and models. Very Easy to Use!
  - \* Just plug in the battery! No need for special adapters, inverters, or voltage matching!
  - \* Powers Mac/PC laptops, internet modems, phones, tablets, lights, or any device < 250W
- **Content:** Battery, wall charger, fast car charger, travel pouch, user manual, storage box.
- **Solar Panel Charger (optional):** Allows travelling in remote areas.
- **Recharge Time:** Less than 12 hours.
- **Unlimited Linking Technology:** No limit on linking up batteries to get more capacity!
- \*\*\*\*\* NOTE: NOT FAA/TSA Compliant for Air Travel (over 160 Wh Limit) \*\*\*\*\*
- **Battery Run Time Chart:** (more details available on Zopec.com)

Heated Humidifier	Heated Tube	CPAP Brand	CPAP Model	Pressure (cm H2O)	Machine Run Time (hours)
OFF	OFF	Philips Respironics	DreamStation	10	64
			DreamStation Go	10	64
		ResMed	AirSense 10	10	50
			AirMini	15	42
On (3)	OFF	Philips Respironics	DreamStation	10	24
				15	18
		ResMed	AirSense 10	10	24
				15	18
On (3)	On (3)	Philips Respironics	DreamStation	10	12
				15	10
		DreamStation BiPap	IPAP = 20	8	
			EPAP = 15		
On (4)	On (71F)	ResMed	AirSense 10	10	12
			15	10	
		AirCurve 10 BiPap	IPAP = 20	8	
			EPAP = 15		