

Introducing the **NEW SMART Travel CPAP Battery**



ZOPEC MEDICAL

Zopec Voyage

dns+ *Wall Outlet on the Go!*



Power Your CPAP Anywhere and Charge Your Phone, Tablet, and Laptop at Lightning Speed!

SMART Travel CPAP Battery. Leave All Tangled Chargers at Home. Just Bring Your USB-C Lightning Wires!

Product Delivery (PD) is the fastest and most versatile SMART charging technology. With PD45W USB-C, it charges 5V/9V/12V/15V/20V devices and at lightning speed, which is 3 - 4 times faster than the regular USB! **SMART PD technology is also built in to Wall Charger, Car Charger, and Solar Chargers!**



Powers All CPAP Brands and Models! No Adapters Needed!

Zopec Voyage has 110V AC Outlet so it works with all CPAP machines. No need to match voltage or figure out the correct adapters!

CPAP Brand	CPAP Model	Pressure (cm H2O)	Heated Humidifier and Heated Tube	Machine Run Time (Hours)
Philips Respironics	DreamStation	10	Both OFF	14
		15	Both OFF	10
	DreamStation Go	10	Both OFF	12
		15	Both OFF	9
ResMed	AirSense 10	10	Both OFF	12
		15	Both OFF	9
	AirMini	10	Both OFF	14
		15	Both OFF	10
Breas HDM	Z1 / Z2	10	Both OFF	16
Somnatics	Transcend	10	Both OFF	14

PRODUCT SPECIFICATIONS

- 2021 New Model (Part#: 20ZDPO2)
- 150 Wh. Long Battery Run Times.
- AC Output: 110V AC/60Hz/100W
- SMART PD45W USB-C Port
- Optional SMART PD60W Solar Charger. Fully charge in 2.5 hours in full sun!
- Optional SMART PD45W Car Charger
- Small Size: 7" x 5" x 1.5" (in)
- Light Weight: 2 lbs.
- Two USB-A Ports: 5V/2.4A
- LED Power Capacity Indicator
- FAA/TSA Compliant for Air Travel
- Recharges less than 4 hours
- Over 500 Charge Cycles
- Product Content: Zopec Voyage Battery, PD39W Wall Charger, Travel Pouch.
- 1 Year Warranty

