

PHILIPS

DreamStation Go
and *DreamStation*

CPAP travel devices
and accessories

For people
on the go



Organized and in style

The Philips Respironics DreamStation and DreamStation Go accessories give you a convenient, compact, discreet and stylish way to travel with your sleep therapy system and help defend against damage including bumps, scratches, dust and more.



DreamStation Go accessories



Heated humidification

Heated humidifiers are designed to add moisture to delivered air with the goal of improving comfort during therapy. DreamStation Go's optional heated humidifier uses tap, bottled or distilled water¹ so travelers don't have the inconvenience of packing or purchasing distilled water when you arrive at your destination.



Overnight battery

Compact and lightweight. DreamStation Go set at 10 cm of pressure displayed an average run time of 13 hours when running off the DreamStation Go's optional overnight battery.

¹In very hard water areas, it is recommended that you use distilled or bottled water.



12 mm micro-flexible tubing

This tubing offers more comfort and greater ease of movement than that of standard tubing.³

³Head-to-head trial of DreamStation 15 mm tubing vs. DreamStation 12 mm micro-flexible tubing, 2016.



DreamStation Go small and medium travel kits

Both kits include a protective case and washable bag for your mask and tube. The medium protective case can accommodate DreamStation Go with either the heated humidifier or battery.

DreamStation accessories



CPAP battery kit¹

Lithium ion battery with integrated uninterruptible power supply (UPS) for DreamStation and System One therapy devices. Perfect for those living on an unreliable power grid or during emergencies.

¹Includes a rechargeable lithium ion battery pack, a CPAP device cable and a battery case.



DreamStation travel case²

Holds your DreamStation³ sleep therapy system mask and tubing in a convenient and compact case.

²DreamStation, mask and tubing not included
³Without humidifier



CPAP travel briefcase⁴

All-in-one carrying option that accommodates most CPAP therapy systems and humidifiers, a laptop and other small items.

⁴DreamStation, power supply and tubing not included



Mask and tubing bedside organizer

The bedside organizer gives you a way to conveniently store your mask and tubing without cluttering a nightstand – both at home and during travel.



DreamStation Go travel CPAP



Mask

Gives you the freedom to choose the mask you want!¹



Air travel

Complies with FAA requirements for in-flight use.⁴



Humidifier

Optional heated humidifier uses tap, bottled or distilled water.²



Charge

DreamStation Go is the only travel CPAP with a USB charging port,³ allowing for easy charging of devices.



Control

No app or mobile device needed to adjust your therapy settings.



Guarantee

Two-year warranty

Whether you're traveling for business or venturing out on the vacation of a lifetime, DreamStation Go provides an easy, reliable and portable CPAP experience for users who refuse to compromise.

¹Masks with standard 22 mm outer diameter connection.

²In very hard water areas, it is recommended that you use distilled or bottled water.

³When compared against HDM ZI, Somnetics Transcend and ResMed's AirMini devices, 2017.

⁴RTCA/DO-160G section 21, category M.

Order your CPAP accessories today

For more information, contact your provider or go to:

sleepapnea.com/dreamstation or sleepapnea.com/dreamstationgo



© 2018 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com/respironics

Caution: U.S. federal law restricts these devices to sale by or on the order of a physician.

RRDPGH MK 12/19/18 MCI 4108569 PN 1140933
1010 Murry Ridge Lane, Murrysville, PA 15668
800 345 6443 · 724 387 4000