

# Interfaces

## FlexiFit™ 432 Full Face Mask



Our second generation Full Face Mask, the **FlexiFit™ 432**, incorporates a soft Foam Cushion that works with the highly contoured Silicone Seal to provide improved sealing technology and a more comfortable fit.

**FlexiFit™** Technology has been incorporated into the contoured Silicone Seal to allow a large range of nasal bridge movement and simple one step auto adjustment fitting. The unique **Glider™** Strap provides freedom of movement while maintaining a seal. The crown strap and Under Chin Design provide extra stability and maintain a fit even with jaw movement. The Easy-Release Strap allows simple removal of the mask and the “Fit and Forget” Headgear ensures no on-going adjustments are required.

The **FlexiFit™ 432** Full Face is supplied fully assembled and is available in Small, Medium or Large.

Product Code	Description	Quantity
HC432	<b>FlexiFit™ 432</b> Full Face Mask includes: <ul style="list-style-type: none"><li>• Silicone Seal (Small, Medium or Large)</li><li>• Foam Cushion (Small, Medium or Large)</li><li>• Headgear with Easy Release Cord</li><li>• Glider™ Strap</li></ul>	1 1 1 1

# FlexiFit™ 432 Mask Spare Parts



## Elbow Kit


Product Code	Description	Quantity
400HC203	FlexiFit™ 431 Elbow & Vented Non Re-Breathing Valve	Each
400HC201	FlexiFit™ 431 Elbow & Non Re-Breathing Valve (Aus only)	Each

**Warning:** Do not use with Non-Invasive Ventilation models.



## Foam Cushion

Product Code	Description	Quantity
400HC004	Foam Cushion FlexiFit™ 432 Small	Each
400HC005	Foam Cushion FlexiFit™ 432 Medium	Each
400HC006	Foam Cushion FlexiFit™ 432 Large	Each




## Foam Cushion & Silicone Seal Combo

Product Code	Description	Quantity
400HC113	Foam & Seal FlexiFit™ 432 Small	Each
400HC114	Foam & Seal FlexiFit™ 432 Medium	Each
400HC115	Foam & Seal FlexiFit™ 432 Large	Each



## Forehead Pads

Product Code	Description	Quantity
900HC423	Forehead Pads - fits FlexiFit™ 405, 407, 431	2 Pack




## Oxygen/Pressure Port Connector

Product Code	Description	Quantity
900HC452	Oxygen/Pressure Port Connector	10 Pack



## Mask Headgear

Product Code	Description	Quantity
400HC309	FlexiFit™ 432 Headgear with Crown Strap	Each



## Mask Kit - No Headgear

Product Code	Description	Quantity
400HC519	FlexiFit™ 432 Small Mask – No Headgear	Each
400HC521	FlexiFit™ 432 Medium Mask – No Headgear	Each
400HC522	FlexiFit™ 432 Large Mask – No Headgear	Each
400HC520	FlexiFit™ 432 Small Mask – No Headgear Australia Only	Each
400HC524	FlexiFit™ 432 Medium - Australia Only	Each
400HC525	FlexiFit™ 432 Large - Australia Only	Each