

## Compliance Letter for In-Flight Use of M-PED

To whom it may concern,


We, Apex Medical Corp., as a manufacturer of Continuous Positive Airway Pressure (CPAP), hereby declare that Apex iCH CPAP Series had been approved by U.S. Food and Drug Administration and Taiwan Food and Drug Administration. The iCH CPAP Series is also in compliance with electromagnetic compatibility in accordance with the current version of RTCA/DO-160, Section 21, Category M and are considered FAA compliant meaning it is safe to be operated on a commercial aircraft.

The U.S. Department of Transportation (DOT) final rule "Nondiscrimination on the Basis of Disability in Air Travel" (effective May 13, 2009) requires that air carriers conducting passenger service must permit someone with a disability, to use a ventilator, respirator, continuous positive airway pressure (CPAP) machine.

Apex iCH CPAP models identified in the table below have been tested and comply with RTCA/DO-160, Section 21, Category M.

iCH Auto	9S-007201, 9S-007200
iCH Prime	9S-007101, 9S-007100
iCH Sense	9S-007001, 9S-007000

Sincerely,



---

Frank Lin

Director of Quality Management Division

April 26, 2016