

Designed **with care**



As the leader in mask technology, ResMed has partnered with pediatric health experts to bring you Pixi™—a pediatric mask designed with children and their families in mind.

A mask designed **specifically for children**

Uniquely constructed using the anthropometric data of children aged two and over, Pixi Pediatric Mask is not merely a scaled down adult mask.

Taking into account the bone structure, facial characteristics and skin sensitivity of children, the mask is designed to minimize pressure on soft faces.

And when combined with ResMed's latest Stellar™ 100 ventilator, you have the ideal pediatric solution that delivers greater comfort and reliability for children over two years of age and weighing more than 13 kg.

Pixi Pediatric Mask is made with both the child and their caregiver in mind.

Every feature is carefully considered to address and improve therapy **comfort, acceptance and overall experience for everyone involved**—so caregivers can be free to focus on emotionally supporting their child on their journey.



Pixi™
PEDIATRIC MASK



Technology Highlight

NEW Headgear designed especially for children*

Three points of adjustment—the mask can be adjusted at the top, side and bottom straps, providing a greater range of adjustment to cater for children of various ages.

Out-of-sight—the headgear is designed to sit away from the child's eyes and ears to reduce irritation and obstruction.

Clear line of sight—not only can the child see clearly while wearing the mask, but carers are also reassured and comforted by being able to see their child's eyes and more of their face.

**Designed using the specific head and facial anthropometric data of children.*

Soft, comfortable cushion

Single-wall, thin silicone cushion minimizes pressure on children's soft faces, without compromising seal

Two tubing positions

Air tubing can be moved to either side of the cushion depending on your child's sleeping position

Spring-flex lightweight tube

stretches and bends easily, so it won't get caught or drag as your child moves around during sleep



Comfortable Breathoprene™ headgear materials

lightly colored for a more reassuring, unobtrusive appearance on the child's face

Side-facing back strap

ergonomically designed from the fitting adult's point of view

Flexible stabilisers

allow the mask to conform to each child's unique facial structure and shape for the most customized fit

Quick-release latch

makes it easy to take the mask off, while retaining desired headgear settings

Emergency latch

Release the latch even faster in an emergency by pulling firmly on the lower headgear strap

Complete System

USA, Asia Pacific,
Australia and New Zealand

Pixi **61030**

Combine the Pixi mask with the latest Stellar 100 device for a system that delivers greater comfort, reliability and effective care.



Stellar 100



Pixi Mask



Pixi™
PEDIATRIC MASK