

Introducing the next generation in Airgel. Our iQ® Blue Nasal Vented and Non-Vented masks feature our patented AIRgel with *Advanced Cushion Technology* and flexible shell with *Custom Fit Technology*™. Together, they provide an optimal fit with maximum comfort. Unlike masks with rigid shells, iQ Blue adjusts to truly fit the contours of your patients' face.

The Smart Choice for a Custom Fit

 **RespsShop**
better sleep, better life.
[Click Here To Buy](#)

Clearest Field of Vision Possible
No forehead pad required.

Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.
Flex, fit, forget.



AIRgel

Maximum comfort, minimum weight. Our patented AIRgel minimizes pressure points and skin irritation.

Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Vented and Non-Vented Versions Available

Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.



What makes our iQ Blue Mask different?

- *Custom Fit Technology* is unique to Sleepnet®. Our soft-shell mask has a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Our *AIR°gel with Advanced Cushion Technology* creates an ergonomic gel cushion that reduces pressure points on the skin.
- Clearest field of vision possible — users can even wear glasses with this mask.
- No forehead pad, yet retains outstanding stability.

About Advanced Cushion Technology (A.C.T.)

Our patented AIR°gel with an ergonomic streamlined design and luxurious feel reduces pressure while maintaining an effective seal.



How does iQ Blue stand apart from the competition?

Features	Sleepnet iQ Mask	ResMed Mirage™ Softgel Mask	Fisher & Paykel Forma Mask	Philips Respironics Comfort Gel Mask
Custom Fit Technology	●			
Advanced Cushion Technology	●			
Minimal forehead contact	●			
Easy universal sizing	●	●	●	●
Breathable headgear	●			
Patented AIR°gel	●			
100% made in the USA	●			

iQ Blue Mask Specifications:

Product Number	Mask Style	Headgear Style	Mask Size	Tube Length	Mask Weight	Mask Dimensions
50647	Vented	3-Point	OS	13"	63.8 g	3.04"D x 2.95"W x 2.60"H
50655	Vented	StableFit	OS	13"	63.8 g	3.04"D x 2.95"W x 2.60"H
50250	Vented	3-Point	OS	6"	63.8 g	3.04"D x 2.95"W x 2.60"H
50255	Vented	StableFit	OS	6"	63.8 g	3.04"D x 2.95"W x 2.60"H
50660	Non-Vented	3-Point	OS	6"	63.9 g	3.08"D x 3.13"W x 2.50"H



Mojo®

Your Mask, Your Way

Our full-face mask also features AIR°gel.



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA



Scan code to view iQ Fitting Tips video.