

Zzz-PAP™

Size, economy, efficiency, integrity...all add up to Zzz-PAP, a breakthrough system that redefines comfort and quality

MINI CPAP UNIT

The Probasics Zzz-PAP™ by PMI delivers performance, simplicity, style, and value in a very small package. The Zzz-PAP offers easy patient setup with a 4-20 cmH₂O pressure range, a comfortable 0-45 minute ramp feature and manual altitude adjustment. PMI also offers the optional Probasics Zzz-PAP Heated Humidifier that quickly and easily integrates to the Zzz-PAP unit providing excellent relief from the symptoms of dryness and congestion associated with CPAP therapy. Weighing less than two pounds and coming standard with an auto-adjusting 100 Volt/240 Volt Universal Power Supply, the compact Zzz-PAP is not only an excellent choice for your bedside but it's perfect for travel as well. Best of all, the Probasics Zzz-PAP is backed by a two-year manufacturer's warranty which we all know is the most necessary feature for providing both the patient and the dealer with peace of mind.



 **Respsshop**
better sleep, better life.
[Click Here To Buy](#)

INTEGRATED HEATED HUMIDIFIER

PMI offers an Integrated Heated Humidifier system designed for use with the Probasics Zzz-PAP™ and Zzz-PAP AUTO units. It relieves nasal irritation and dryness which result from constant airflow by supplying the required humidity to the upper airway. With a reliable and safe design featuring overheat protection and a water chamber safety switch, the Probasics Zzz-PAP Integrated Heated Humidifier can drastically improve the patient's comfort level and compliance



**The Zzz-PAP™ and
Zzz-PAP AUTO units
mount directly on the
Integrated
Heated Humidifier**

Zzz-PAP AUTO™

Simply outstanding!
...the ultimate in simplicity, functionality and value.

AUTO CPAP UNIT

Much like the the Probasics Zzz-PAP™, the Probasics Zzz-PAP™ AUTO delivers performance, simplicity and value in a very small package. The difference is the Zzz-PAP™ AUTO also offers an Auto CPAP mode with automatic altitude adjustment as well as compliance and efficacy monitoring with download to a Mini-SD Card. The Zzz-PAP™ AUTO's flow sensing technology utilizes a special honeycomb structure to create a constant flow resistance for minimizing pressure swing. The net result is a quiet, smooth and effective unit that provides effective Auto-CPAP therapy. The Probasics Zzz-PAP™ AUTO also offers features such as automatic shut-off and a clock with a functional alarm. Each unit comes standard with a mini SD card with SD adaptor (USB-SD card reader sold separately), USB cord, and Zzz-Ware Compliance and Efficacy Software.



Mini SD Card & SD Adaptor USB-SD Card Reader

Zzz-WARE Compliance & Efficacy Software: Captures, Displays, Graphs and Prints Total Operating Time of the Unit, Compliance Data and Event Data (AI, HI, SI).

ORDERING INFORMATION

- 7520 Zzz-PAP Core Package**
(Unit Complete with Heated Humidifier)
- 7521 Zzz-PAP Unit Only**
- 7620 Integrated Heated Humidifier only**
- 7601 Replacement Humidifier Chamber**
- 7510 Replacement Filters (Pack of 10)**

- 7540 Zzz-PAP AUTO Core Package**
(Unit Complete with Heated Humidifier)
- 7541 Zzz-PAP AUTO (Unit only)**
- 7545 USB - SD Card Reader**

Note: All Zzz-PAP and Zzz-PAP AUTO units come complete with Operator's Manual, Carrying Bag, 6' Tube and AC Power Cord.

SPECIFICATIONS

Zzz-PAP™

C-PAP and Auto C-PAP Units

	Zzz-PAP	Zzz-PAP AUTO
Pressure Range	4.0 - 20.0 cmH ₂ O (0.5 cmH ₂ O increment)	4.0 - 20.0 cmH ₂ O (0.5 cmH ₂ O increment)
Ramp Time	0 - 45 min (5 minute increments)	0 - 45 min (5 minute increments)
Altitude Compensation	Manual Adjustment	Auto Adjustment
Leak Compensation	N/A	Yes
Compliance Monitoring	N/A	Yes
Events Detection (AI/HI/SI)	N/A	Yes
Download Capacity	N/A	Yes
Clock w/Alarm Function	N/A	Yes
Device Operating Meter	9999.9 hr	9999.9 hr
Heated Humidifier	Optional	Optional
Sound Level	30 dBA @ 10 cmH ₂ O	30 dBA @ 10 cmH ₂ O
Dimension (W x D x H)	5.8"x5"x3.9"	5.8"x5"x3.9"
Weight	1.77 lbs.	1.77 lbs.
Power Supply	100V-240V, 50/60Hz, Universal Power Supply	100V-240V, 50/60Hz, Universal Power Supply
Operating Temperature	41°F to 95°F	41°F to 95°F
Environment Humidity	15% to 95% RH	15% to 95% RH
Warranty	Two years	Two years

Zzz-PAP™

Integrated Heated Humidifier

- Integrates with Probasics **Zzz-PAP** and **Zzz-PAP AUTO** units
- Easy to use rotary knob for heating level adjustment
- Dishwasher safe water chamber has an easy to remove base plate for cleaning
- Safe & reliable with overheat protection and water chamber safety switch
- Indicators clearly show when **Power is ON** and **Heater is ON**

Power Supply	100V-240V, 50/60Hz (Universal)
Water Capacity	450 ml
Heater Settings	1 - 6 (Heater plate temperature 86°F to 140°F)
Warm-Up Time	< 30 minutes
Pressure Drop	0.2 cmH ₂ O @ 60LPM
Humidity Output	10 - 40 mgH ₂ O/L
Indicator	Power On indicator, Heater On indicator
Dimension (W x D x H)	Heater: 6.8" x 10" x 1.6" Water Chamber: 4.5" x 5.9" x 4.2"
Weight	Heater: 0.84 lb Water Chamber: 0.53 lb
Warranty	One year

PROBASICS™
by PMI

Probasics™

15 South Main Street

Second Floor, Unit F

Marlboro, NJ 07746

Tel. 800.438.0503

Fax 732.683.1945

www.pmiusa.biz

All specifications are subject to change without notice.

Economies of scale

Zzz-PAP AUTO™

Zzz-PAP™

PROBASICS™
by PMI

